

Single-Use Cannulas

Our new range of sterile, single-use cannulas is specially crafted to enhance safety and simplicity in liposuction and infiltration procedures. Designed for ease of use without the need for complex sterilization processes, our cannulas prioritize enhanced patient safety. With a variety of shapes available, they serve the needs of current liposuction techniques with precision.

ADVANTAGES

- + minimize the risk of infections and cross-contamination and improve the patients safety.
- + ready to use – the time consuming process for cleaning, disinfection and sterilization is obsolete and increases the efficiency.
- + include a variety of shapes and fulfill the requirements of the end-user.

+ NOUVAG AG


St. Gallerstrasse 25
9403 Goldach
Switzerland
T +41 (0) 71 846 66 00
info@nouvag.com

+ NOUVAG GmbH

Schulthaißstrasse 15
78462 Konstanz
Germany
T +49 (0) 7531 1290-0
F +49 (0) 7531 1290-12
info-de@nouvag.com

nouvag.com

©NOUVAG AG • 67018EN • V20240820 • All rights reserved



NOUVAG⁺

Single-Use Cannulas

Aspiration Cannulas



SINGLE PORT			
REF		4643	4644
Cannula Ø	mm	2.0	3.0
Length	mm	150	300



THREE PORT			
REF		4647	4648 4649
Cannula Ø	mm	2.0	3.0 4.0
Length	mm	200	300 300



MERCEDES			
REF		4653	4654 4655 4656
Cannula Ø	mm	2.0	3.0 3.0 4.0
Length	mm	150	250 300 250



TRIPPLE MERCEDES	
REF	4659
Cannula Ø	mm 5.0
Length	mm 300



DOUBLE BASKET			
REF		4662	4663 4664
Cannula Ø	mm	4.0	5.0 5.0
Length	mm	300	300 250



SATTLER					
REF		4666	4667	4668	4669 4670
Cannula Ø	mm	2.0	3.0	3.0	3.0 4.0
Length	mm	150	200	300	250 300
Holes		30	30	30	30 30



BAYONETTE		
REF		4673 4674
Cannula Ø	mm	3.0 4.0
Length	mm	300 300
Holes		10 10



ONE HOLE SPATULA	
REF	4675
Cannula Ø	mm 2.0
Length	mm 150



30 HOLES				
REF		4676	4677	4678 4679
Cannula Ø	mm	2.0	2.0	2.0 3.0
Length	mm	150	250	300 250
Holes		30	30	30 30



TWO PORT			
REF		4645	4646
Cannula Ø	mm	2.5	3.0
Length	mm	200	300




ONE CENTRAL, TWO LATERAL			
REF		4650	4651 4652
Cannula Ø	mm	3.0	4.0 5.0
Length	mm	300	300 300



DOUBLE MERCEDES		
REF		4657 4658
Cannula Ø	mm	4.0 5.0
Length	mm	250 300



SINGLE BASKET		
REF		4660 4661
Cannula Ø	mm	3.0 4.0
Length	mm	250 300



CANDY CANE	
REF	4665
Cannula Ø	mm 4.0
Length	mm 300



TOLEDO V		
REF		4671 4672
	mm	2.0 3.0
Length	mm	150 250

Infiltration Cannulas

THREAD CONNECTION
internal thread

